



**DIREKTORAT PENGELOLAAN LABORATORIUM,
FASILITAS RISET DAN KAWASAN SAINS DAN
TEKNOLOGI**

Kawasan Sains dan Teknologi Soekarno, Jalan Raya Bogor Km. 46
Cibinong, Jawa Barat 16911

FORMULIR DETAIL PENGUJIAN

Nama :
Instansi :
*) S1/S2/S3/Industri/Umum/Internal BRIN
Alamat :
No Hp/Telp :
E-mail :

**CHECKLIST PARAMETER
ANALISIS OLIGOSAKARIDA (HPLC)**

No	Parameter	
1	Jumlah Sample	:
2	Konsentrasi Sample	:
3	Fase Gerak (Mobile Phase)	:
	Detector	: <u>Refractive Index Detector (RID) – 20A</u>
4	Kolom	: Aminex HPX-42A, 300 x 7.8 mm
5	Kolom Oven Temperature	:
6	Kolom Flow Mode	:
7	Komponen Target	: <input type="checkbox"/> Xylooligosakarida (X1-X6) <input type="checkbox"/> Mannooligosakarida (M1-M6) <input type="checkbox"/> Maltooligosakarida (M1-M6)
9	Catatan Lain-lain	: